



## **Biological Air Sampler TRBS-601**

[www.labdistributor.com](http://www.labdistributor.com) | [info@labdistributor.com](mailto:info@labdistributor.com)

## Description

Biological Air Sampler TRBS-601 is highly accurate and has a fast sampling rate. It has a sample selection amount setting that can be adjustable from 0.01 to 2.0 liters, with a measuring range of 10 to 6000 L. It comes with an LCD that shows the sampling quantity 3.53 CFM/min, time, and other parameters. It has a collection efficiency of 95%. It is programmable and simple to operate. To reduce puffball superposition and animal culturing error, the sampling hole has many small holes.

## Features:

The sampling hole has a multitude of different holes to reduce puffball quantum state and animal culture count error  
It has sampling quantities ranging from 0.01 to 2.0 liters  
The LCDs show the sampling quantity, sampling time, and other parameters  
Saving of sampling quantity, sampling time, and other parameters in pages, with a maximum storage capacity  
Sculpting that is simple to use  
It has a simple to exchange substrate utensil, simply take down sampling hole to exchange substrate utensil

## Specifications :

<b>Collection Efficiency</b>	0.95
<b>Power Supply</b>	AC & DC, chargeable battery: DC7.4V, continuous work for 4h after charged
<b>Sampling Quantity</b>	3.53CFM/min
<b>Sampling Range</b>	10 to 6000L
<b>Speed Crash</b>	16m/s
<b>Speed Of Sampling Hole</b>	0.4m/s, basically the same as that of cleanroom
<b>Volume</b>	?120×300 mm
<b>Weight</b>	2.5Kg



Main Street Kansas City MO 64108, USA  
Email: [info@labdistributor.com](mailto:info@labdistributor.com) | Website: [www.labdistributor.com](http://www.labdistributor.com)