



## **Fume Hood TRFH-601**

[www.labdistributor.com](http://www.labdistributor.com) | [info@labdistributor.com](mailto:info@labdistributor.com)

## Description

Fume Hood TRFH-601 is a ventilation cabinet with a solvent disinfectant storage tank made from steel. It is secure and dependable. Internal air expellees at a rate of 0.5 m/s, resulting in a high ventilation effectiveness and a ventilation capacity of 450 to 1600 m<sup>3</sup>/h. It is having with a 1.5 M fume hood. It has a liquid crystal controller for the display. The upper section is fitted with an exhaust fan system and a moving glass door.

## Features:

Water, electricity, gas, and air conditioning are all interrogatively into the cabinet  
Inside, there is a multifunctional socket  
A stainless steel water tank with liquid antiseptic is installation on the bottom  
It is secure and dependable  
The ventilation capacity is 450 to 1600 m<sup>3</sup>/h, and the air circulation effectiveness is high  
It has a sliding glass door that can be open and closed

## Specifications :

<b>Dimension Of Working Area (W1×D1×H1)</b>	1330×700×600mm
<b>Discharging Speed</b>	within 0.3 m to 0.6m/s adjustable
<b>Inlet Pressure</b>	>0.5Pa
<b>Overall Dimension W×D×H</b>	1500×760×1950mm
<b>Power</b>	800W
<b>Power Supply</b>	AC, 220 V/50Hz
<b>Weight</b>	180 Kg



Main Street Kansas City MO 64108, USA  
Email: [info@labdistributor.com](mailto:info@labdistributor.com) | Website: [www.labdistributor.com](http://www.labdistributor.com)