



Digital Heating Water Bath TRHB-602

www.labdistributor.com | info@labdistributor.com

Description

Digital Heating Water Bath TRHB-602 is a microprocessor controller with an inner digital display showing a temperature range of (RT+5) -300 °. It holds with a capacity of 6 L. It is designed with a stainless steel interior with an aluminum ring cover.

Features:

It is equipped with a microprocessor control
mechanism with a stainless steel interior
Designed with a digital display inner temperature
It is covered with an aluminum ring
It is a thermostatic device

Specifications :

Capacity	6 L
Dimension	37x23x24 cm
Power	600 W
Temperature Fluctuation	±0.5 °
Temperature Homogeneity	±1 °
Temperature Range	(RT+5) -300 °
Voltage	110V 60Hz /220V 50 Hz
Weight	4.5 kg
Working Size	30x16x12 cm



Main Street Kansas City MO 64108, USA
Email: info@labdistributor.com | Website: www.labdistributor.com