



Digital Heating Water Bath TRHB-604

www.labdistributor.com | info@labdistributor.com

Description

Digital Heating Water Bath TRHB-604 is a microprocessor controller with an inner digital display showing a temperature range of (RT+5) -300 °. It holds with a capacity of 17 L. It is designed with a stainless steel interior with an aluminum ring cover.

Features:

It is equipped with a microprocessor control mechanism with a stainless steel interior
Designed with a digital display inner temperature
It is covered with an aluminum ring
It is a thermostatic device

Specifications :

Capacity	17 L
Dimension	53x37x24 cm
Power	1500 W
Temperature Fluctuation	±1.5 ?
Temperature Homogeneity	±1 ?
Temperature Range	(RT+5) -300 ?
Voltage	110V 60Hz /220V 50 Hz
Weight	7.5 kg
Working Size	46x30x12 cm



Main Street Kansas City MO 64108, USA
Email: info@labdistributor.com | Website: www.labdistributor.com