



Portable Ozone Generator TROG-602

www.labdistributor.com | info@labdistributor.com

Description

Portable Ozone Generator TROG-602 is a lightweight and portable device. It can use for a variety of sterilization and purification treatments. It applies to both air and water, and it operates on the tube type-high voltage discharged gap principle. Designed with stainless steel and therefore can start producing O₃ at a rate of 5 g/h. It performs and has a long service life.

Features:

It is having with good performance and long life

It is equipped with a stainless steel chamber

This applies to the use of air and water due to the operating principle of tube-type high voltage discharge paths

It has a lightweight and portable device

It can generate O₃ ranges

It is having with an ozone concentration

Intended with a high voltage ionization

It is both industry and household type

It is an automatic control time type device

It is having with a cool wind of Byo

Specifications :

Application	Both can use air and water
O3 range	5h/g
Overall size	500 ×160× 220 mm
Packing Size	550×220×270 mm
Power	220V/ 50HZ



Main Street Kansas City MO 64108, USA
Email: info@labdistributor.com | Website: www.labdistributor.com