



Water Steam Distiller TRWD-604

www.labdistributor.com | info@labdistributor.com

Description

Water Steam Distiller TRWD-604 is a stainless steel water tower system that is corrosion resistant and resistant to aging. Distillation of generating equipment, condenser, evaporator, electric heating pipe, and control unit. It has an automatic, intelligent control system. It is equipped with a control panel and has a button-type mode. Designed with a water level window, which can clearly show the potresidual water level.

Features:

It is a water tower mechanism made from steel
It is corrosion-resistant and age resistant
It has an automatic and intelligent control system
designed with a water level window that displays the residual water level
The heating component is made up of immersion heating tubes that are long-lasting

Specifications :

Output	5 L/ hour
Package	540×270×490mm
Power	AC220 V/ 4.5 KW
Weight	9 Kg



Main Street Kansas City MO 64108, USA
Email: info@labdistributor.com | Website: www.labdistributor.com